mega **CRYL HOT + DON**

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For further Informations contact us !

mega CRYL HOT + DON

hot curing resin

mega CRYL HOT + DON

is a hot curing resin on a Methylmethacrylatbasis for the manufacture of total and partial prostheses, done in a pressure-process. Good flowability, high mechanical values, secure Polymerisation

Working instructions: The plaster surface must very clean from wax, cover the plaster surface 2x with **measol** separator. Put Liquid into the mixing-bowl, then add the powder and mix in a fast motion, until there is no more powder absorbed. After about 10 minutes rising time the paste can be kneaded for approx 45 minutes (by 18-24°C)

Mixing ratio: 10 parts powder : 4 parts liquid (or as desired)

Polymerisation: Put flask in cold or max. 30-40°C warm water, heat up water to boiling (98-100°C) make cooking along 25-30 minutes, and cool down in the water.

Short time-polymerisation: Normal (regular) pieces into hot (80°C) water, turn off heat and wait 15-20 minutes. Then let water come to a boil again and let the Polymerisation take place for 10-15 min.

Colours: pink, pink-veined, light-pink, light-pink-veined, pink-opak, clear
Packing: 0,5 kg / 1 kg / 5 kg / 20 kg # 0,5 Lit. / 1 Lit. / 5 Lit.

	Lot number		manufaturer
	Expiry date		Xi, irritant
	Read instructions for use		F, Highly inflammable
	Recommended storage temperature		

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

megadental is certified in accordance to ISO 9001 and guarantees efficient quality through maintaining a deliberate quality management system. Product properties as mentioned in this instruction can only be guaranteed when the procedural steps are carried out properly. The user is responsible for correct handling of the product. Any defective results cannot be accounted for, because megadental has no influence on the subsequent processing by the user. Any claims for damages can solely be granted for the value of our goods.

according to medical device directive 93/42 EEC annex IX class **Ila** for fixed and removable dentures
according to DIN EN ISO 20795 Type 1 Group 1

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