


SPLIT CAST




Installation of Split Cast Kit




Split Cast Kit Réf QM 2000




2 screws , 2 magnets
2 split cast plates
and 2 retention disks






1 : remove the nut of the screw










2 : Fixe the screws on articulator




3 : *Fixe the magnets on articulator*









Réf QM 2000

Split Cast Kit

2 screws , 2 magnets
2 split cast plates
and 2 retention disks



Réf QM 2010

**Split Cast Plates (x 24)
with
Retention Disks (x 24)**



Réf QM 2020

Retention Disks (x 50)

Prices Ex-works factory

April 1st 2016

		Unit Prices in Euros		
		1 to 4 units (same réf.)	5 to 12 units (same réf.)	+ 12 units (same réf.)
<u>SPLIT CAST SYSTEM</u>				
QM2000	2 screws, 2 magnets, 2 split cast plates and 2 retention disks	35,29 €	32,08 €	28,87 €
QM2010	Split Cast plates (X24) retention disks (X24)	39,79 €	36,17 €	32,55 €
QM2020	Retention Disks (X50)	18,61 €	16,92 €	15,23 €