




fag dentaire

CLUSES





COMPARATIVE SUPPORT





fag dentaire



Quick Master LAB



Quick Master B2 M

Type	Arcon	Arcon
Sagittal Condylar Inclination	25° or 40°	0° to 70°
Bennett Angle	0° - 10° - 15° - 20°	0° - 10° - 15° - 20°
Other Settings	I.S.S	I.S.S - Propulsion
Use With Facebow	Yes	Yes
Price	141 €	216 €



Asa 5000



Asa 5010



Asa 5020



Asa 5030

Type	No - Arcon	No - Arcon	No - Arcon	No - Arcon
Sagittal Condylar Inclination	30°	30°	30°	30°
Bennett Angle	7°	7°	7°	7°
Other Settings	No	No	No	No
Use With Facebow	No	No	No	Yes
Price	~ 71 €	~ 89 €	~ 94 €	~ 107 €



PROTHE Plus evo 3



PROTHE Plus evo 5



PROTHE Plus evo 7



PROTHE Plus evo 9

Type	Arcon	Arcon	Arcon	Arcon
Sagittal Condylar Inclination	30°	- 25° to 75°	- 25° to 75°	- 25° to 75°
Bennett Angle	15°	15°	0° to 30°	0° to 30°
Other Settings	No	No	ISS , retrusion , protrusion , distraction	ISS , retrusion , protrusion , distraction
Use With Facebow	Yes	Yes	Yes	Yes
Price	~ 608 €	~ 1240 €	~ 1520 €	~ 1870 €



AMANNGIRRBACH



ARTEX BN



ARTEX CN



ARTEX CT

Type	No - Arcon	No - Arcon	No - Arcon
Sagittal Condylar Inclination	35°	35°	- 15° to 60°
Bennett Angle	15°	0° to 20°	0° to 20°
Other Settings	No	No	No
Use With Facebow	No	No	Yes
Price	~ 286 €	~ 394 €	~ 685 €



AMANNGIRRBACH



ARTEX CP



ARTEX CR

Type	Arcon	Arcon
Sagittal Condylar Inclination	- 20° to 60°	- 20° to 60°
Bennett Angle	-5° to 30°	-5° to 30°
Other Settings	Distraction	ISS , retrusion , protrusion , distraction
Use With Facebow	Yes	Yes
Price	~ 982 €	~ 1374 €



Stratos 100



Stratos 300

Type	Arcon	Arcon
Sagittal Condylar Inclination	30°	0° to 60°
Bennett Angle	15° and 30°	0° to 30°
Other Settings	No	ISS - Protrusion
Use With Facebow	Yes	Yes
Price	~ 422 €	~ 1044 €



SAM[®]



SAM SE



SAM 2P



SAM 2PX



SAM 3

Type	Arcon	Arcon	Arcon	Arcon
Sagittal Condylar Inclination	0° to 65°	0° to 65°	0° to 65°	0° to 65°
Bennett Angle	?	0° to 45°	0° to 45°	0° to 45°
Other Settings	No	ISS , retrusion , protrusion	ISS , retrusion , protrusion	ISS , retrusion , protrusion
Use With Facebow	Yes	Yes	Yes	Yes
Price	~ 410 €	~ 804 €	~ 985 €	~ 1394 €



Type	No - Arcon	No - Arcon	Arcon	Arcon
Sagittal Condylar Inclination	?	- 10° to 70°	30°	0° to 60°
Bennett Angle	?	0° to 25°	15°	0° to 30°
Other Settings	No	No	No	No
Use With Facebow	No	Yes	Yes	Yes
Price	~ 85 €	~ 212 €	~ 141 €	~ 259 €



Baumann

Dental GmbH



ArTO 1



ArTO 2



ArTO 3

<i>Type</i>	<i>No - Arcon</i>	<i>No - Arcon</i>	<i>No - Arcon</i>
<i>Sagittal Condylar Inclination</i>	<i>35°</i>	<i>35°</i>	<i>- 15° to 60°</i>
<i>Bennett Angle</i>	<i>12°</i>	<i>0° to 20°</i>	<i>0° to 20°</i>
<i>Other Settings</i>	<i>No</i>	<i>No</i>	<i>No</i>
<i>Use With Facebow</i>	<i>Yes</i>	<i>Yes</i>	<i>Yes</i>
<i>Price</i>	<i>~ 298 €</i>	<i>~ 358 €</i>	<i>~ 478 €</i>



mälzer dental

PRÄZISIONSTEILE FÜR DIE ZAHNTECHNIK



Mal 12100 - 12200



Mal 12110 - 12210



Mal 12300



Mal 12310

Type	No - Arcon	No - Arcon	No - Arcon	No - Arcon
Sagittal Condylar Inclination	35°	- 15° to 60°	35°	-15° to 60°
Bennett Angle	0° to 25°	0° to 25°	0° to 25°	0° to 25°
Other Settings	No	No	No	No
Use With Facebow	Yes	Yes	Yes	Yes
Price	~ 289 €	~ 389 €	~ 319 €	~ 419 €



Automark



Mark 310



Anamark Plus

<i>Type</i>	<i>Arcon</i>	<i>Arcon</i>	<i>Arcon</i>
<i>Sagittal Condylar Inclination</i>	<i>25°</i>	<i>25°</i>	<i>- 0° to 60°</i>
<i>Bennett Angle</i>	<i>7°</i>	<i>15°</i>	<i>7°</i>
<i>Other Settings</i>	<i>No</i>	<i>No</i>	<i>ISS</i>
<i>Use With Facebow</i>	<i>Yes</i>	<i>Yes</i>	<i>Yes</i>
<i>Price</i>	<i>~ 370 €</i>	<i>~ 501 €</i>	<i>~ 613 €</i>



Denar



Combi II




Mark 320



Mark II

Type	Arcon	Arcon	Arcon
Sagittal Condylar Inclination	0° to 60°	0° to 60°	0° to 90°
Bennett Angle	15°	15°	5° to 15°
Other Settings	No	No	ISS
Use With Facebow	Yes	Yes	Yes
Price	~ 716 €	~ 726 €	~ 783 €



Mark 330



Track II



D5A

Type	<i>Arcon</i>	<i>Arcon</i>	<i>Arcon</i>
Sagittal Condylar Inclination	<i>0° to 60°</i>	<i>0° to 90°</i>	<i>0° to 60°</i>
Bennett Angle	<i>0° to 15°</i>	<i>3° to 30°</i>	<i>0° to 30°</i>
Other Settings	<i>ISS</i>	<i>ISS - Protrusion</i>	<i>ISS</i>
Use With Facebow	<i>Yes</i>	<i>Yes</i>	<i>Yes</i>
Price	<i>~ 875 €</i>	<i>~ 901 €</i>	<i>~ 1831 €</i>